


# CrystalQuick 96 Well Plate, 2 µl round drop well, 3 drop well, LBR, treated, Greiner 609820 - 10/cs

Sitting drop vapor diffusion crystallization



### **Product codes:**

Reference: HR 3090

## **Product gallery:**





### **Product description:**

- SBS specifications for high throughput crystallization and are fully compatible with robotic equipment
- 96 well sitting drop plate
- 3 drops/well
- Round drop wells
- 0.1 to 1.9 µL drop volume
- Standard profile (14.4 mm plate height)
- 320 µL maximum reagent volume
- Cyclic polyolefin low birefringence (LBR) material for polarized crystal imaging, exceptional transparency, high chemical resistance, and low water absorption

The plates can be sealed using Crystal Clear Sealing Film (HR3-609), 3 inch wide Crystal Clear Sealing Tape (HR4-506) or ClearSeal Film (HR4-521).

To quickly compare features and select the Greiner CrystalQuick plate that best fits your needs please refer to the Greiner CrystalQuick Plate Summary pdf on this page.



Per maggiori informazioni visita il sito https://hamptonresearch.com/

# **Product features:**

CRF - TIPO: CrystalQuick 96 Well