



24% w/v Polyethylene glycol 3,500 - 24% w/v solution - 200 ml - 200 ml

Crystallization screen for intact monoclonal antibodies, monoclonal antibody fragments, & proteins less soluble at low ionic strength

Low Ionic Strength Screen™

Product codes:

Reference: HR 2519

Product gallery:



Product description:

- SBS microplate footprint & standard well spacing
- Optically clear wells provide for superior crystal imaging
- Well layout makes crystal harvesting easier
- 3 wells are ideal for optimizing protein concentration, additive screens, drop ratio and combinatorial experiments
- Rounded reservoir corners to prevent reservoir (MPD) creep
- Plate Material: Cyclic Olefin
- UV transmissible and low birefringence

Art Robbins Instruments has developed the ideal family of plates for crystallography applications: the Intelli-Plate. All versions of the Intelli-Plate have the SBS microplate footprint and standard well spacing and are ideal for manual or automated processing.

The Intelli-Plate 96-3 LVR (Art Robbins 102-0001-03 Hampton Research HR3-183 and HR3-185) is a Low Volume Reservoir, 3 drop well version of the Intelli-Plate 96-2 Original (Art Robbins 102-0011-00 Hampton Research HR3-297 and HR3-299).

Reservoir volume: 100 µl max.



Typical reservoir volume 50 μ l.
Drop well volume: 1 μ l max.
Plates are made of cyclic olefin.
UV transmissible and low birefringence.

The plates can be sealed using Crystal Clear Sealing Film (HR3-609), 3 inch wide Crystal Clear Sealing Tape (HR4-506) or ClearSeal Film (HR4-521).

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: Low Ionic Strength Screen