



18 mm Mounted CryoLoop Sampler - 20 micron (5 of each loop diameter size 0.1 - 1.0 mm) - 30/pk

Cryocrystallography

Product codes:

Reference: HR 4918

Product gallery:



Product description:

- Ready to use 18 mm length MicroTube with CryoLoop
- Synthetic CryoLoop on stainless steel MicroTube
- Complete range of loop diameters from 0.025 - 1.0 mm
- Compatible with CrystalCap, CrystalCap ALS, CrystalCap ALS HT and CrystalCap Spine HT

18 mm Mounted CryoLoops with 20 micron diameter nylon. Nylon loops are staked to hollow, stainless steel MicroTubes that are used to mount, freeze, and secure the crystal during cryocrystallographic procedures and X-ray data collection. The MicroTube on the 18 mm Mounted CryoLoops is 18 mm in length and ready to insert into CrystalCap Magnetic and CrystalCap HT. The 20 micron diameter nylon loops show minimal diffraction, are thin for fast cryo cooling of the crystal, strong, and aerodynamic.

For use with CrystalCap, CrystalCap ALS, CrystalCap ALS HT and CrystalCap SPINE HT systems. Not compatible with the CrystalCap Copper, CrystalCap Copper HT, or CrystalCap ALS HT.



Shop name

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: 18 mm Mounted CryoLoop - 20 micron