




Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted

Cf 25pz - Piastra sub-skirted 96x0.2ml (RP) Spec. ABI/Life T. - frosted. Cf 25pz

Product codes:

Reference: 12 893 0681



Product gallery:



Product description:

Piastra da 0.2ml semi skirted. Utilizzabile su sistemi PCR e PCR Real Time a profilo regolare. Questa piastra è inoltre perfettamente compatibile con sistemi ABI e Life Technologies

- Colore bianco
 - Regular Profile
 - Sub skirted
 - Miscela in PP specifico per applicazioni in PCR
 - DNA, RNA, DNase, RNase e Pyrogen FREE
- Extra Low Binding
- Autoclavabili

Product features:

PSC - TYPE: plates PCR
PSC - VOLUME: 0.2ml
PSC - FORM: sub-skirted
PSC - OPENING: normal
PSC - WELLS: 96
PSC - BOTTOM: a V
PSC - CLOSING: mat
PSC - CLOSING: via caps



PSC - CLOSING: through film
PSC - APPLICATIONS: PCR Real-Time
PSC - APPLICATIONS: PCR
PSC - APPLICATIONS: molecular biology
PSC - AUTOCLAVABLE: you
PSC - COLOUR: frosted
PSC - PURITY: DNA / RNA / DNase / RNase / Pyrogen FREE
PSC - STERILITY: no
PSC - MATERIAL: PP
PSC - COMPLIANCE: not subject to CE
PSC - TEMPERATURE MAX: 100°C
PSC - MINIMUM TEMPERATURE: 30°C
PSC - PACKAGE: bulk
PSC - MAIN FEATURES: piastra PCR 96 pozzetti da 0.2 specifiche ABI e Stratagene per PCR e PCR Real Time
semi-skirted - FROSTED
PSC - OPTIONAL ACCESSORIES: mats sealants
PSC - OPTIONAL ACCESSORIES: film sealants
PSC - OPTIONAL ACCESSORIES: the strip of caps
PSC - SUITABLE FOR: PCR
PSC - SUITABLE FOR: PCR Real-Time
PSC - SUITABLE FOR: storage
PSC - SUITABLE FOR: sample preparation
PSC - PROFILE: adjust
PSC - COMPATIBLE WITH: Stratagene Mx3005P™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 6 Flex 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 7 Flex 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems SimpliAmp™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Solid3 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ 96-Well Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems Veriti™ Dx 96-Well Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ViiA™ 7 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Stratagene Mx3000P™ 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems MiniAmp Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7300 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 2700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7500 Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBLE WITH: Applied Biosystems ABI PRISM 7000 SDS 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ABI PRISM 7700 SDS 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 5700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems MiniAmp Plus Thermal Cycler 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 3 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 310 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3100 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3130 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 3730 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems 7900HT Fast Real-Time PCR System 96 x 0.2ml Standard
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 2720 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Dual 96-Well 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Aluminum 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9700 Gold 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems ProFlex™ PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 12k Flex 96 x 0.2ml



PSC - COMPATIBLE WITH: Applied Biosystems QuantStudio™ 5 Real-Time PCR System 96 x 0.2ml
PSC - COMPATIBLE WITH: Applied Biosystems GeneAmp PCR System 9600 96 x 0.2ml
PSC - CAPACITY: 0.2ml
PSC - WALLS: thin 0.35 mm
PSC - TREATMENT SIGNAL: walls of frosted
PSC - SIGNAL-STRENGTH: 125
PSC - INTERFERENCE SIGNAL: 70
PSC - RESISTANT: freezing at -20°C
PSC - EVAPORATION SAMPLE: the lesser of 3%
PSC - BINDING: DNA minor 0.3%
PSC - BINDING: proteins less 0.8%
PSC - GRIGLIA ALFANUMERICA: si