





CrystalQuick 96 Well Plate, 4 μ L square drop well, 1 drop well, untreated, low profile, Greiner 609171 - 20/cs

Sitting drop vapor diffusion crystallization

Product codes:

Reference: HR 3302

Product gallery:



Product description:

- SBS specifications for high throughput crystallization and are fully compatible with robotic equipment
- 96 well sitting drop plate
- 1 drop/well
- Square drop wells
- 0.1 to 4.1 μ L drop volume
- Low profile (8.0 mm plate height)
- 140 μ L maximum reagent volume
- Cyclic polyolefin low birefringence (LBR) material for polarized crystal imaging, exceptional transparency, high chemical resistance, and low water absorption

The plates can be sealed using Crystal Clear Sealing Film (HR3-609), 3 inch wide Crystal Clear Sealing Tape (HR4-506) or ClearSeal Film (HR4-521).



To quickly compare features and select the Greiner CrystalQuick plate that best fits your needs please refer to the Greiner CrystalQuick Plate Summary pdf on this page.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: CrystalQuick 96 Well