



Natrix 2 06 / Natrix HT E6, 185 ml

0.05 M Sodium cacodylate trihydrate pH 6.0, 35% v/v Tacsimate pH 6.0, 0.001 M Spermine,

, Individual Natrix - Natrix 2 and Natrix HT screen reagents for screening optimization and production



Product codes:

Reference: HR 291706

Product gallery:



Product description:

185 ml volume – Sterile filtered

– Formulated in Type 1+ ultrapure water: 18.2 megaohm-cm resistivity at 25°C, 5 ppb Total Organic Carbon, bacteria free (1 Bacteria (CFU/ml)), pyrogen free (0.03 Endotoxin (EU/ml)), RNase-free (0.01 ng/mL) and DNase-free (4 pg/μL)

HR2-918-**

** Refers to the reagent number in the kit.

For example, Natrix reagent number 1 = Natrix HT reagent A1 = HR2-918-01

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: Individual Natrix - Natrix 2 - Natrix HT Reagents