




Grid Screen PEG/LiCl - 10 ml. tube format - 10 ml

Primary or secondary, PEG & Lithium chloride versus pH grid crystallization screen for biological macromolecules

Product codes:

Reference: HR 2217

Product gallery:



Grid Screen™
PEG/LiCl



Product description:

- SBS microplate footprint & standard well spacing
- Optically clear wells provide for superior crystal imaging
- Well layout makes crystal harvesting easier
- 2 wells with small and large drop well
- Rounded reservoir corners to prevent reservoir (MPD) creep
- Plate Material: Cyclic Olefin
- UV transmissible and low birefringence

Art Robbins Instruments has developed the ideal family of plates for crystallography applications: the Intelli-Plate. All versions of the Intelli-Plate have the SBS microplate footprint and standard well spacing and are ideal for manual or automated processing.

The Intelli-Plate 96-2 LVR (Art Robbins 102-0001-00 Hampton Research HR3-143 and HR3-145) is a Low Volume Reservoir version of the Intelli-Plate 96-2 Original (Art Robbins 102-0011-00 Hampton Research HR3-297 and HR3-299).

Reservoir volume: 100 µl max.

Drop well volume: 1 µl max.



Plates are made of cyclic olefin.
UV transmissible and low birefringence.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

Product features:

CRF - TIPO: Grid Screen PEG/LiCl