



# Brodo bile verde brillante al 2 % (BRILA) BAC ISO-4831, ISO-4832, 500 g

Brodo bile verde brillante al 2 % (BRILA) BAC ISO-4831, ISO-4832, 500 g



## Codes produits :

Référence 71 899 1156

## Galerie de produits :

## Description du produit :

Brodo Bile Verde Brillante 2 % (BBVB) BAC ISO-4831, ISO-4832  
TARIC: 13023100

Homogeneous, free flowing and Beige-greenish.

If any changes physically, discard product.

Appearance Brilliant green

pH (25°C) 7.0 – 7.4pH

Concentration of use 40 g/L

Preparation:

Dissolve 40 gr in 1 L of water with heat and boil for 1 min until complete dissolution.

Sterilize at 121°C for 15 minutes

(AVOID OVERHEATING)

Once opened avoid hydration



Storage temperature 2 °C – 25 °C

**caractéristiques du produit :**

RGT - GRADO: BAC - Bacteriology Grade

RGT - FAMIGLIA: Lifekem

TAR - TYPE: réactifs

TAR - APPLICATIONS EN CHIMIE: bottiglia di plastica

RGT - CAPACITA': 500gr