



Green Gel Loading Buffer with DNA Stain - 5 x 1,8ml

Green Gel Loading Buffer with DNA Stain

Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain



Codes produits :

Référence JN PCR-275

Galerie de produits :

Description du produit :

Green Gel Loading Buffer with DNA Stain - Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain

5 x 1,8ml

Shipping Conditions: shipped on gel packs

Storage Conditions: store at -20 °C
store dark

Shelf Life: 12 months

eClass Nummer: 32160502

Product Web Link: [Product Web Link](#)

Product Data Sheet: [Download Data Sheet](#)

Product Safety Data Sheet: [Download Product Safety Data Sheet](#)

caractéristiques du produit :

TAR - TYPE: Reagenti

TAR - TYPE: Polymerase chain reaction PCR and reverse transcriptase polymerase chain reaction RT PCR products

TAR - DURÉE D'EXPIRATION: 12 months

TAR - TRANSPORT: shipped on gel packs



RGT - CAPACITA': 5 x 1

RGT - CONSERVAZIONE: Store at -20 °C

TAR - NUMÉRO CAS: -