



# Orange Gel Loading Buffer with DNA Stain - 5 x 1,8ml

Orange Gel Loading Buffer with DNA Stain

Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain



**Codes produits :**

Référence JN PCR-276

**Galerie de produits :**

## **Description du produit :**

Orange Gel Loading Buffer with DNA Stain - Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain

5 x 1,8ml

**Shipping Conditions:** shipped on gel packs

**Storage Conditions:** store at -20 °C  
store dark

**Shelf Life:** 12 months

**eClass Nummer:** 32160502

Product Web Link: [Product Web Link](#)

Product Data Sheet: [Download Data Sheet](#)

Product Safety Data Sheet: [Download Product Safety Data Sheet](#)

## **caractéristiques du produit :**

TAR - TYPE: Reagenti

TAR - TYPE: Polymerase chain reaction PCR and reverse transcriptase polymerase chain reaction RT PCR products

TAR - DURÉE D'EXPIRATION: 12 months

TAR - TRANSPORT: shipped on gel packs



---

RGT - CAPACITA': 5 x 1

RGT - CONSERVAZIONE: Store at -20 °C

TAR - NUMÉRO CAS: -