



# PEGRx 2 29 / PEGRx HT G5, 185 ml

6% v/v Tacsimate pH 6.0, 0.1 M MES monohydrate pH 6.0, 25% w/v Polyethylene glycol 4,001,  
, Individual PEGRx 1, PEGRx 2, and PEGRx HT screen reagents for screening, optimization and production

# PEGR<sub>x</sub><sup>TM</sup>

## Codes produits :

Référence HR 298429

## Galerie de produits :

## Description du produit :

- 185 ml volume
- Sterile filtered
- Formulated in Type 1+ ultrapure water: 18.2 megaohm-cm resistivity at 25°C, 5 ppb Total Organic Carbon, bacteria free (1 Bacteria (CFU/ml)), pyrogen free (0.03 Endotoxin (EU/ml)), RNase-free ( 0.01 ng/mL) and DNase-free ( 4 pg/μL)

HR2-982-\*\*

\*\* = Refers to the reagent number in the kit.

For example, PEGRx 1 reagent number 1 = PEGRx HT reagent A1 = HR2-982-01

For example, PEGRx 2 reagent number 1 = PEGRx HT reagent E1 = HR2-984-01

For PEGRx 2 reagent 48(H12) 3% w/v Dextran sulfate sodium salt Mr 5,000, 0.1 M BICINE pH 8.5, 15% w/v Polyethylene glycol 20,000, please order HR2-984-48A for PEGRx 2 kits lot 208404, 208405, 208406, and order HR2-428-51B for PEGRx 2 kits lot 208401, 208402, 208403 and 2084zz where zz >6.

For PEGRx HT reagent H12(48) 3% w/v Dextran sulfate sodium salt Mr 5,000, 0.1 M BICINE pH 8.5, 15% w/v Polyethylene glycol 20,000, please order HR2-428-51A for PEGRx HT kits lot



---

208609-04-04, 208610-04-04, 208611-05-04, 208612-05-05, 208613-06-06, 208614-06-06, 208615-07-06. For all other PEGRx HT kit lots order HR2-428-51B.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

**caractéristiques du produit :**

CRF - TIPO: Individual PEGRx 1 - PEGRx 2 - PEGRx HT Reagents