

Brodo Muller Kauffmann con novobiocina e verde brillante BAC ISO-6579, 500 g

Brodo Muller Kauffmann con novobiocina e verde brillante BAC ISO-6579, 500 g

**Product codes:**

Reference: 71 899 3166

Product gallery:**Product description:**

Brodo Muller Kauffmann con novobiocina e verde brillante BAC ISO-6579
TARIC: 13023100

Homogeneous, free flowing and white.

If any changes physically, discard product.

Appearance Green-blue

pH (25°C) 7.8 – 8.2pH

Concentration of use 89.53 g/L

Preparation:

Dissolve 89.53 gr on 1 L with heat and boil for 1 min until complete dissolution.

Cool to 45-50°C and aseptically add 20 ml of IODI-R0A POIO-00A solution (20g IODI 25g POIO in 100ml of sterile distilled water)

(DO NOT AUTOCLAVE)

(AVOID OVERHEATING)

Once opened avoid hydration

Storage temperature 2 °C – 25 °C

Product features:

RGT - GRADO: BAC - Bacteriology Grade

RGT - FAMIGLIA: Lifekem

RGT - TIPO: Chimici

RGT - PACKAGING: bottiglia di plastica

RGT - CAPACITA': 500gr