

Orange Gel Loading Buffer with DNA Stain - 5 x 1,8ml

Orange Gel Loading Buffer with DNA Stain

Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain



Product codes:

Reference: JN PCR-276

Product gallery:

Product description:

Orange Gel Loading Buffer with DNA Stain - Loading buffer for agarose or polyacrylamide gels with green-fluorescent DNA stain

5 x 1,8ml

Shipping Conditions: shipped on gel packs

Storage Conditions: store at -20 °C
store dark

Shelf Life: 12 months

eClass Nummer: 32160502

Product Web Link: [Product Web Link](#)

Product Data Sheet: [Download Data Sheet](#)

Product Safety Data Sheet: [Download Product Safety Data Sheet](#)

Product features:

RGT - TIPO: Reagenti

RGT - TIPO: Polymerase chain reaction PCR and reverse transcriptase polymerase chain reaction RT PCR products

RGT - DURATA SCADENZA: 12 months

RGT - TRASPORTO: shipped on gel packs

RGT - CAPACITA': 5 x 1

RGT - CONSERVAZIONE: Store at -20 °C

RGT - CAS NUMBER: -