


# CrystalQuick 96 Well Plate, 4 $\mu$ l square drop well, 3 drop well, untreated, Greiner 609101 - 10/cs

Sitting drop vapor diffusion crystallization

## Product codes:

Reference: HR 3192

## Product gallery:



## Product description:

- SBS specifications for high throughput crystallization and are fully compatible with robotic equipment
- 96 well sitting drop plate
- 3 drops/well
- Square drop wells
- 0.1 to 4.1  $\mu$ L drop volume
- Standard profile (14.4 mm plate height)
- 320  $\mu$ L maximum reagent volume
- Polystyrene material

The plates can be sealed using Crystal Clear Sealing Film (HR3-609), 3 inch wide Crystal Clear Sealing Tape (HR4-506) or ClearSeal Film (HR4-521).

To quickly compare features and select the Greiner CrystalQuick plate that best fits your needs

please refer to the Greiner CrystalQuick Plate Summary pdf on this page.

Per maggiori informazioni visita il sito <https://hamptonresearch.com/>

**Product features:**

CRF - TIPO: CrystalQuick 96 Well